



CRYO VIALS

High Performance Lab Polyester (B-461)
FreezerBondz™ (B-490)
Nylon Cloth (B-499)


Part No.	Vial Size	Diagram	Reference	Labels per Cartridge	Dimensions – mm			Install Style
					A	B	C	
LONG-TAIL STYLE – text covered for solvent resistance								
60941	.6 ml	D	XSL-73-461	225	9.53	45.72	19.05	Long-Tail
60940	.6 ml	D	XSL-72-461	225	12.70	44.45	19.05	Long-Tail
60942	1.0 to 8.0 ml	D	XSL-74-461	175	12.70	55.88	19.05	Long-Tail
60944	1.0 to 8.0 ml (larger text area)	D	XSL-86-461	150	15.24	66.68	25.40	Long-Tail
60943	1.0 to 8.0 ml (larger text area)	D	XSL-75-461	150	25.40	66.68	25.40	Long-Tail
60933	12 to 14 ml	D	XSL-119-461	100	38.10	95.25	31.75	Long-Tail
SHORT-TAIL STYLE – lowest profile, easier handling								
60935	.6 ml	D	XSL-123-461	500	9.53	28.58	19.05	Short-Tail
60934	.6 ml	D	XSL-122-461	500	12.70	28.58	19.05	Short-Tail
60940	1.0 to 8.0 ml	D	XSL-72-461	225	12.70	44.45	19.05	Short-Tail
60938	1.0 to 8.0 ml (larger text area)	D	XSL-126-461	225	15.24	45.72	25.40	Short-Tail
60937	1.0 to 8.0 ml (larger text area)	D	XSL-125-461	225	25.40	44.45	25.40	Short-Tail
60939	12 to 24 ml	D	XSL-129-461	175	38.10	57.15	31.75	Short-Tail
CONTINUOUS WRAP STYLE – use w/vial templates in LABXPERT, no window to see liquid level								
60930	.6 ml	⊙	XC-375-461	30 ft roll	9.1 m	9.53		Short Overlap
60931	.6 to 2.0 ml	⊙	XC-500-461	30 ft roll	9.1 m	12.70		Short Overlap
60927	1.0 to 8.0 ml	⊙	XC-1000-461	30 ft roll	9.1 m	25.40		Short Overlap
60928	4.0 to 8.0 ml	⊙	XC-1250-461	30 ft roll	9.1 m	31.75		Short Overlap
FREEZERBONDZ™ MATERIAL								
60903	.6 ml	▲	X-112-490	500	31.75	9.53		Short Overlap
60909	1.0 to 8.0 ml	▲	X-126-490	225	45.72	15.24		Short Overlap
60908	1.0 to 8.0 ml	★	X-125-490	225	25.40	44.45		Short Overlap
60905	14 to 15 ml	★	X-113-490	175	38.10	55.88		Short Overlap
NYLON CLOTH MATERIAL								
60911	.6 ml	▲	X-130-499	350	20.96	9.53		Free Standing
60950	.6 ml	⚡	X-118-499	350	25.40	9.53	9.53 dia	Free Standing
60949	.6 ml	⚡	X-120-499	350	25.40	12.70	9.53 dia	Free Standing
60913	1.0 to 8.0 ml	▲	X-131-499	350	25.40	12.70		Free Standing
60948	1.0 to 8.0 ml	⚡	X-117-499	350	25.40	12.70	11.18 dia	Free Standing
60915	12 to 24 ml	★	X-132-499	225	31.75	38.10		Free Standing



VIAL TOP LABELS

Nylon Cloth (B-499)

Part No.	Vial Size	Diagram	Reference	Labels per Cartridge	Dimensions – mm			Color
					A	B	C	
60465	.6 to 1.5 ml	✘	X-82-499	350	9.530			BK on WT
60923	.6 to 1.5 ml	✘	X-82-499-YL-BK	350	9.530			BK on YL
60921	.6 to 1.5 ml	✘	X-82-499-OR-BK	350	9.530			BK on OR
60922	.6 to 1.5 ml	✘	X-82-499-RD-BK	350	9.530			BK on RD
60950	.6 to 1.5 ml	⚡	X-118-499	350	25.40	9.530	9.530 dia	BK on WT
60949	.6 to 1.5 ml	⚡	X-120-499	350	25.40	12.70	9.530 dia	BK on WT



IDENTILAB™ LABORATORY LABELLING SOFTWARE

Create files, templates and more with the easy-to-use wizards in Brady's IdentiLab™ Laboratory Labeling Software. The streamlined interface and hundreds of parts and graphics categorised by common lab applications makes creating legible, durable and accurate labels quick and easy.

Compatible with a variety of Brady printers including LABXPERT™ Laboratory Labeling System, TLS 2200® Thermal Labeling System, TLS-PC Link™ Thermal Labeling System, all Bradyprinter™ models, and all Windows-based printers.

- Easy-to-use wizards to create new files and templates
- Superscript and subscript characters
- Database import capabilities with refresh feature to ensure data integrity
- Automatic text formatting
- 2-D bar code support Multi-label data entry view
- Die cut and continuous label supplies

Part No.	Reference	Cartridge
51961	IDL	Identilab V1.0

